


Electrophysiology stimuli for DIGAP

Supplementary Information

Supplementary document for Koos et al., Automatic deep learning driven label-free image guided patch clamp system for human and rodent in vitro slice physiology manuscript.

Stimulation waveform to test electrode resistance:



Stimulation waveform to analyze the physiological properties of patched neurons:

Name & Description	Length	Amplitude	Repeat interval
IVCC - Incremented long pulse	800 ms	-100 pA	0 s

The figure displays a diagram of an incremented long pulse (IVCC) showing a series of horizontal steps. A scale bar indicates 20 pA (vertical) and 100 ms (horizontal).

increment: 20 pA

Batch commands for controlling the HEKA EPC amplifiers:

PROTOCOL "RESET"

This protocol is used for resetting the oscilloscope window in PatchMaster and switch the amplifier to voltage clamp mode.

```
Command ( 0.000s): " O DispTrace      0; Trace 1"
Command ( 0.000s): " O ResetY"
Command ( 0.000s): " O ResetX"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O DispTrace      1; Trace 2"
Command ( 0.000s): " O ResetY"
Command ( 0.000s): " O ResetX"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O DispTrace      0; Trace 1"
Command ( 0.000s): " E Ampl3          TRUE"
Wait     ( 0.000s): abs 100.0ms
Command ( 0.000s): " E Ampl1          FALSE"
Wait     ( 0.000s): abs 100.0ms
Command ( 0.000s): " E Mode           3; Whole Cell"
Wait     ( 0.000s): abs 100.0ms
Switch   ( 0.000s): "Oscilloscope"
Switch   ( 0.000s): "Amplifier"
Wait     ( 0.000s): abs 100.0ms
```

PROTOCOL "bHunt"

This protocol is used for creating a new group in the file structure in PatchMaster, set the test pulse, set the gain of the amplifier, calculate the offset potential of the electrode, set the filter, and call the "shortpulse" protocol. This protocol is used when the pipette is in the bath. The "shortpulse" protocol is running during the cell hunting phase, and the recorded current signals measured with the NI board to calculate the pipette tip resistance.

```
Command ( 0.000s): " @ File           "New Group""
Chain    ( 0.000s): "RESET", return
Command ( 0.000s): " E Reset"
Command ( 0.000s): " E Mode           3; Whole Cell"
Command ( 0.000s): " E PulseAmp      -5.0mV"
```

```

Command ( 0.000s): " E PulseDur 5.0ms"
Command ( 0.000s): " E CSlow 100.00pF"
Command ( 0.000s): " E RSeries 20.0MOhm"
Command ( 0.000s): " E Gain 7; 0.5 mV/pA -> medium
range"
Command ( 0.000s): " E AutoZero"
Command ( 0.000s): " E Filter2 7.4kHz"
Command ( 0.000s): " E PulseOn TRUE"
Command ( 0.000s): " E SaveRpip"
Series ( 0.000s): "shortpulse","",""
Wait ( 0.000s): abs 100.0ms
Command ( 0.000s): " E Gain 10; 1.0 mV/pA ->
medium range"
Command ( 0.000s): " E Ampl2 TRUE"
Wait ( 0.000s): abs 100.0ms
Command ( 0.000s): " E Ampl1 TRUE"
Wait ( 0.000s): abs 100.0ms
Switch ( 0.000s): "Oscilloscope"
Switch ( 0.000s): "Amplifier"
Wait ( 0.000s): abs 100.0ms
Command ( 0.000s): " N Store TRUE"

```

PROTOCOL "bBreakin"

This protocol is used for compensating the fast capacitive transient of the pipette. The protocol is called if the pipette tip resistance higher than 1 GOhm ie. When the gigaseal is formed.

```

Chain ( 0.000s): "RESET", return
Command ( 0.000s): " E Mode 3; Whole Cell"
Command ( 0.000s): " E Gain 10; 5 mV/pA"
Command ( 0.000s): " E AutoCFast"
Command ( 0.000s): " E AutoCFast"
Switch ( 0.000s): "Amplifier"

```

PROTOCOL "aBreakin"

##This protocol is used for the cell capacitance compensation and for Rs compensation. It is called when the whole cell configuration is achieved.

```

Chain ( 0.000s): "RESET", return
Command ( 0.000s): " E Mode 3; Whole Cell"
Command ( 0.000s): " E Gain 11; 5 mV/pA"
Command ( 0.000s): " E CSlow 10.00pF"
Command ( 0.000s): " E RSeries 20.0MOhm"
Command ( 0.000s): " E AutoCSlow"
Command ( 0.000s): " E AutoCSlow"
Command ( 0.000s): " E Gain 11; 10 mV/pA"
Switch ( 0.000s): "Amplifier"
Chain ( 0.000s): "IVCC"

```

PROTOCOL "IVCC"

This protocol is used for switching the amplifier to current clamp mode and start the recording. From this point, all the electrophysiological recordings are automatized in PatchMaster.

```
Command ( 0.000s): " N Store TRUE"
Command ( 0.000s): " E Ampl1 FALSE"
IF ( 0.000s): AD-7 > -100.00p
    Command ( 0.000s): " E Mode 4; C-Clamp"
    Command ( 0.000s): " E IHold 0.0pA"
END_IF
Command ( 0.000s): "E Mode 4; C-Clamp"
Command ( 0.000s): " O DispTrace 0; Trace 1"
Command ( 0.000s): " O ResetY"
Command ( 0.000s): " O ResetX"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O YScaleInc"
Command ( 0.000s): " O YScaleInc"
Wait ( 0.000s): abs 150.0ms
Value ( 0.000s): Value-4 = 0.0000
Value ( 0.000s): Value-3 = 0.0000
SetPgf ( 0.000s): PgfParam-1 = 0.0000
REPEAT ( 0.000s): sweeps 0.000s
    Sweep ( 0.000s): "Long square","", ""
    IF ( 0.000s): Value-3 = 0.0000
        IF ( 0.000s): Value-1 > 1.0000
            Value ( 0.000s): Value-4 = Value-2
            Value ( 0.000s): Value-3 = 1.0000
        END_IF
    END_IF
    IF ( 0.000s): Value-1 > 16.000
        BREAK ( 0.000s): repeat
    END_IF
END_REPEAT
Value ( 0.000s): Value-4 MUL 1.0000G
SetPgf ( 0.000s): PgfParam-1 = Value-4
Chain ( 0.000s): "RHEOBASE"
Command ( 0.000s): " E Mode 3; Whole Cell"
Command ( 0.000s): " E VHold -70.0mV"
Switch ( 0.000s): "Amplifier"
```

See more: https://www.heka.com/support/tutorials/tutorials_down/pm_tutorial.pdf (Chapter 14).